

Instruction Manual to application of program for quantitative detection of *Legionella pneumophila* DNA using AmpliSens® *Legionella pneumophila*-FRT PCR kit

Quantitative characterized calibrators are used for quantification of *Legionella pneumophila* DNA and Internal Control STI-338 DNA copies in a test tube.

To perform the calculation one should take account of Internal Control DNA losses in order to calculate the real concentration of *Legionella pneumophila* DNA in each test sample of water.

To perform the calculation one should take account of **water concentration level**. **Therefore the water pretreatment should be strictly performed in accordance with the *Instruction Manual* to the reagent kit.**

Calculation of concentration values of *Legionella pneumophila* DNA ($C_{DNA Lp}$) in 1 l of water is performed in accordance with the following formula:

$$C_{DNA Lp} \text{ (copies/l)} = Q_{DNA Lp} / Q_{IC STI-338} * C_{IC STI-338} * 2, \text{ where:}$$

$C_{DNA Lp}$ (copies/l) is the quantity of *Legionella pneumophila* DNA copies in 1 l of water sample,

$Q_{DNA Lp}$ (copies/ml) is the calculated quantity of *Legionella pneumophila* DNA copies in 1 ml of a test sample,

$Q_{IC STI-338}$ (copies/ml) is the calculated quantity of Internal Control STI-338 DNA copies in 1 ml of the Internal Control in a test sample,

$C_{IC STI-338}$ (copies/ml) is the number copies of Internal Control STI-338 DNA in 1 ml of Internal Control (specified in the *Important Product Information Bulletin*,

2 is the recalculation coefficient which take account of the volume changes during filtration.

For quantitation each water sample should be tested in two runs
NOTE: (beginning with DNA extraction stage). In such case the result is given as an average value of two obtained values.

The *Program for Calculation of Legionella DNA Quantity.xls* can be used to calculate the concentration of *Legionella pneumophila* DNA.

Procedure

1. Open *Program for Calculation of Legionella DNA Quantity.xls*
2. Save file with **Book Microsoft Excel** extension, specify the date of analysis.
3. In the **Quant. results – Cycling A.FAM / Cycling A. Green** window of Rotor-Gene program copy **Name, Type, Ct, Given Conc (copies/ml), Calc Conc (copies/ml)**

columns and paste them into similar columns of **FAM / Green** field (**Excel** program).

4. In the **Quant. results – Cycling A.JOE / Cycling A. Yellow** window of **Rotor-Gene** program copy **Name, Type, Ct, Given Conc (copies/ml), Calc Conc (copies/ml)** columns and paste them into similar columns of **JOE / Yellow** field (**Excel** program).
5. Enter the concentration of value Internal Control STI-338 specified in the *Important Product Information Bulletin* into all rows of **C Internal Control STI-338** column (**Excel** program).
6. **Concentrations of *Legionella pneumophila* DNA (copies/l of tested water sample) calculated automatically in accordance with the given formula will be displayed in the C DNA Lp column of the program.**